# BOND RNAscope® Brown Detection

Catalog No: DS9815

Leica Biosystems Newcastle Ltd  
Balliol Business Park West  
Benton Lane  
Newcastle Upon Tyne NE12 8EW  
United Kingdom  
+44 191 215 4242

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## Catalog No: DS9815

### Intended Use

This detection system is for *in vitro* diagnostic use.

The BOND RNAscope® Brown Detection enables the user to perform chromogenic *in situ* hybridization (CISH) with Advanced Cell Diagnostics (ACD) proprietary RNA ISH probes. Performance claims for compatible probes have not been established. Probes used in conjunction with BOND RNAscope® Brown Detection should be validated by the user in accordance with local laws and regulations. The BOND RNAscope® Brown Detection is intended for staining sections of formalin-fixed, paraffin-embedded (FFPE) tissue on the Leica BOND-III system.

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls. They should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

The BOND RNAscope® Brown Detection must be used with laboratory best practice in the use of tissue controls. For assurance, laboratories should stain each test sample in conjunction with positive, negative and other tissue specific controls as needed. Appropriate control tissue should be run on the same slide as test samples to ensure optimum quality control.

### Summary and Explanation

*In situ* hybridization techniques can be used to demonstrate the presence of RNA in tissue and cells (see "Using BOND Reagents" in your BOND user documentation).

The BOND RNAscope® Brown Detection utilizes a novel branched DNA-like technology to amplify target probe hybridization signals. This detection system allows visualization of target RNA molecules through chromogenic conversion of DAB by HRP, resulting in a brown chromogenic stain. The BOND RNAscope® Brown Detection enables chromogenic RNA ISH to be performed on the automated BOND-III system. Using the BOND RNAscope® Brown Detection in combination with the BOND-III automated system reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

### Reagents Provided

The reagents provided are sufficient for 12 individual staining runs on the BOND-III, or a maximum of 60 slides.

To achieve a maximum of 60 slides from this detection system, slides must be batched in quantities of 5 or greater per slide staining assembly. Batching in quantities of less than 5 will result in fewer stained slides.

RNAscope® Rinse (27 mL)

RNAscope® Rinse (27 mL)

Hematoxylin (9 mL)

DAB Part 1 (1.2 mL)

DAB Part B (22 mL)

RNAscope® Bluing (9 mL)

RNAscope® H<sub>2</sub>O<sub>2</sub> (9 mL)

RNAscope® AMP 1 DAB (18 mL)

RNAscope® AMP 2 DAB (18 mL)

RNAscope® AMP 3 DAB (18 mL)

RNAscope® AMP 4 DAB (18 mL)

RNAscope® AMP 5 DAB (18 mL)

RNAscope® AMP 6 DAB (18 mL)

RNAscope® Protease (12 mL)

### Dilution and Mixing

BOND RNAscope® Brown Detection is optimized for use on the BOND-III system. Reconstitution, mixing, dilution, or titration of these reagents is not required.

### Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and *in situ* hybridization staining using the BOND-III system.

### Storage and Stability

Store at 2–8 °C. Do not freeze. Do not use after the expiration date indicated on the tray handle label. Return to 2–8 °C immediately after use.

There are no obvious signs to indicate instability of this product, therefore positive and negative controls should be run simultaneously with unknown specimens (refer to "Quality Control" in the "Using BOND Reagents" section of your BOND user documentation).

If unexpected staining is observed that cannot be explained by variations in laboratory procedures, and a problem with the detection system is suspected, contact your local distributor or the regional office of Leica Biosystems immediately.

Storage conditions other than those specified above must be verified by the user<sup>1</sup>.

## Precautions

- This detection system is intended for *in vitro* diagnostic use.

<b>DAB Part 1</b> Contains Ethylene Glycol (>90%) and 66 mM (<10%) 3,3' diaminobenzidine tetrahydrochloride hydrate. GHS07: Exclamation mark. GHS08: Health hazard. Signal words: Danger.	H302: Harmful if swallowed. H341: Suspected of causing genetic defects. H350: May cause cancer.	P201: Obtain special instructions before use. P202: Do not handle until all safety precautions have been read and understood. P280: Wear protective gloves/protective clothing/eye protection/face protection. P264: Wash hands thoroughly after handling. P270: Do not eat, drink or smoke when using this product. P301+312: IF SWALLOWED: Call a POISON CENTER or doctor if you feel unwell. P330: Rinse mouth. P308+313: IF exposed or concerned: Get medical advice/attention. Restricted to professional users.
<b>RNAscope® AMP 1 DAB</b> Contains Formamide (<30%). GHS08: Health hazard. Signal words: Danger.	H360D: May damage the unborn child.	P201: Obtain special instructions before use. P202: Do not handle until all safety precautions have been read and understood. P280: Wear protective gloves/protective clothing/eye protection/face protection. P308+313: IF exposed or concerned: Get medical attention. P501: Dispose of contents/container to hazardous or special waste collection point. Restricted to professional users.
<b>RNAscope® AMP 2 DAB</b> Contains a mixture of 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one (3:1) (<0.0025%). GHS07: Exclamation mark. Signal words: Warning.	H317: May cause an allergic skin reaction.	P261: Avoid breathing mist. P272: Contaminated work clothing should not be allowed out of the workplace. P280: Wear protective gloves/protective clothing/eye protection. P302+352: IF ON SKIN: Wash with plenty of water. P333+313: If skin irritation or rash occurs: Get medical advice/attention. P362+364: Take off contaminated clothing and wash it before reuse. P501: Dispose of contents/container to hazardous or special waste collection point.
<b>RNAscope® AMP 3 DAB</b> Contains Formamide (<30%). GHS08: Health hazard. Signal words: Danger.	H360D: May damage the unborn child.	P201: Obtain special instructions before use. P202: Do not handle until all safety precautions have been read and understood. P280: Wear protective gloves/protective clothing/eye protection/face protection. P308+313: IF exposed or concerned: Get medical attention. P501: Dispose of contents/container to hazardous or special waste collection point. Restricted to professional users.
<b>RNAscope® AMP 4 DAB</b> Contains a mixture of 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one (3:1) (<0.0025%). GHS07: Exclamation mark. Signal words: Warning.	H317: May cause an allergic skin reaction.	P261: Avoid breathing mist. P272: Contaminated work clothing should not be allowed out of the workplace. P280: Wear protective gloves/protective clothing/eye protection. P302+352: IF ON SKIN: Wash with plenty of water. P333+313: If skin irritation or rash occurs: Get medical advice/attention. P362+364: Take off contaminated clothing and wash it before reuse. P501: Dispose of contents/container to hazardous or special waste collection point.
<b>RNAscope® AMP 5 DAB</b> Contains a mixture of 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one (3:1) (<0.0025%). GHS07: Exclamation mark. Signal words: Warning.	H317: May cause an allergic skin reaction.	P261: Avoid breathing mist. P272: Contaminated work clothing should not be allowed out of the workplace. P280: Wear protective gloves/protective clothing/eye protection. P302+352: IF ON SKIN: Wash with plenty of water. P333+313: If skin irritation or rash occurs: Get medical advice/attention. P362+364: Take off contaminated clothing and wash it before reuse. P501: Dispose of contents/container to hazardous or special waste collection point.
<b>RNAscope® AMP 6 DAB</b> Contains a mixture of 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one (3:1) (<0.0025%). GHS07: Exclamation mark. Signal words: Warning.	H317: May cause an allergic skin reaction.	P261: Avoid breathing mist. P272: Contaminated work clothing should not be allowed out of the workplace. P280: Wear protective gloves/protective clothing/eye protection. P302+352: IF ON SKIN: Wash with plenty of water. P333+313: If skin irritation or rash occurs: Get medical advice/attention. P362+364: Take off contaminated clothing and wash it before reuse. P501: Dispose of contents/container to hazardous or special waste collection point.

<b>RNAscope® Protease</b>	H317: May cause an allergic skin reaction.	P261: Avoid breathing mist. P272: Contaminated work clothing should not be allowed out of the workplace. P280: Wear protective gloves/protective clothing/eye protection. P302+352: IF ON SKIN: Wash with plenty of water. P333+313: If skin irritation or rash occurs: Get medical advice/attention. P362+364: Take off contaminated clothing and wash it before reuse. P501: Dispose of contents/container to hazardous or special waste collection point. Restricted to professional users.
Contains a mixture of 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one (3:1) (<0.0025%). GHS07: Exclamation mark. Signal words: Warning.		

- To obtain a copy of the Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems Web site, [www.LeicaBiosystems.com](http://www.LeicaBiosystems.com)
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions<sup>2</sup>. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user<sup>1</sup>.
- Do not mix reagents from different detection systems.

### Instructions for Use

BOND RNAscope® Brown Detection is developed for use on the automated BOND-III system with BOND ancillary reagents and user selected RNA ISH probes. The default staining protocol for the BOND RNAscope® Brown Detection on BOND is \*RNAscope DAB ISH Protocol B. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances.

Operating parameters for application of the detection system reagents on the BOND Processing Module have been optimized by Leica Biosystems. These can be displayed by following the instructions in your BOND user documentation.

### Product Specific Limitations

The Detection System is designed for use with proprietary probes; laboratories should contact ACD for access to compatible RNA probes. Probes used in conjunction with BOND RNAscope® Brown Detection should be validated by the user in accordance with local laws and regulations.

The BOND RNAscope® Brown Detection was developed for use on the BOND-III automated system using \*RNAscope DAB ISH Protocol B. Operating parameters for application of the detection system reagents on the BOND Processing Module have been optimized by Leica Biosystems. These can be displayed by following the instructions in your BOND user documentation. BOND RNAscope® Brown Detection has been optimized at Leica Biosystems for use with BOND ancillary reagents. Users who deviate from recommended test procedures must appropriately validate the system and accept responsibility for interpretation of patient results under these circumstances.

Performance conditions may vary due to variation in tissue fixation and the effectiveness of RNA enhancement, and must be determined empirically. Positive and negative reagent controls should be used when optimizing retrieval conditions.

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls. They should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

The BOND RNAscope® Brown Detection must be used with laboratory best practice in the use of tissue controls. For assurance, laboratories should stain each test sample in conjunction with positive, negative and other tissue specific controls as needed.

Appropriate control tissue should be run on the same slide as test samples to ensure optimum quality control.

Contact your local distributor or regional office of Leica Biosystems for further information.

### Troubleshooting

Refer to reference 3 for remedial action.

If the test result does not correspond to the expected results with the use of controls, the test should be repeated.

If the staining result is not as expected, and you wish to troubleshoot performance of the instrument and detection system independently, your local Leica representative can provide specific protocols. The detection kit must be used in conformance to the package instructions and within the shelf life indicated on the product itself.

### Further Information

Further information on *in situ* hybridization with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

### Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order Code M29-P.
3. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In situ Hybridization. A practical approach.* 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

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### Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order Code M29-P.
3. Wilkinson DG. The theory and practice of in situ hybridization. In: Wilkinson DG. (ed.) In situ Hybridization. A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

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## Bibliografi

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order Code M29-P.
3. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In situ Hybridization. A practical approach*. 2nd Edition. New York: Oxford University Press, 1998, pp.18-20.

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## Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order Code M29-P.
3. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In situ Hybridization. A practical approach*. 2nd Edition. New York: Oxford University Press, 1998. pp.18–20.

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## Literatúra

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order Code M29-P.
3. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In situ Hybridization. A practical approach*. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

## Dátum vydania

22 apríl 2020

Leica Biosystems Newcastle Ltd  
Balliol Business Park West  
Benton Lane  
Newcastle Upon Tyne NE12 8EW  
United Kingdom  
+44 191 215 4242

Leica Biosystems Canada  
71 Four Valley Drive  
Concord, Ontario L4K 4V8  
Canada  
+1 800 248 0123

Leica Biosystems Inc  
1700 Leider Lane  
Buffalo Grove IL 60089  
USA  
+1 800 248 0123

Leica Biosystems Melbourne Pty Ltd  
495 Blackburn Road  
Mt Waverley VIC 3149  
Australia  
+61 2 8870 3500