



ANTIBODIES

- Exclusively manufactured by BioVendor
- Unique antibodies
- Validated for common applications - ELISA, WB, IHC, etc.
- High specificity antibodies
- Affinity purified
- Primary antibodies - monoclonal, polyclonal
- Antibody conjugates
- Custom antibody development
- High lot-to-lot consistency
- Lyophilized in PBS
- No additives affecting further conjugation

BIOVENDOR - MANUFACTURING
PRIMARY ANTIBODIES SINCE 1995

Please find more information on our 200+ antibodies on www.biovendor.com/antibodies



MOLECULE	ANTIBODY	APPLICATION*	NOTE
Adiponectin Human	Goat polyclonal	WB, ELISA	affinity purified
	Rabbit polyclonal	WB, ELISA, IHC	affinity purified
	Sheep polyclonal	WB, ELISA, IHC	affinity purified
	Mouse monoclonal, Clone: 5E4	WB, ELISA, IHC	from supernatant
FABP4 Human	Goat polyclonal	WB, ELISA, IHC	affinity purified
	Rabbit polyclonal	WB, IHC	affinity purified
	Sheep polyclonal	ELISA	affinity purified
Cystatin C Human	Rabbit polyclonal	ELISA	affinity purified
	Sheep polyclonal	WB, ELISA, IHC	affinity purified
Club Cell protein human	Rabbit polyclonal	WB, ELISA, IHC	affinity purified
COMP Human	Mouse monoclonal, Clone: 16F12	WB, ELISA, IHC	from serum free supernatant
	Mouse monoclonal, Clone: 17C10	ELISA, IHC	from serum free supernatant
	Mouse monoclonal, Clone: 12C4	WB, IHC	from serum free supernatant
Fetuin A Human	Goat polyclonal	WB, ELISA, IHC, Immunoturbidimetry	affinity purified
	Rabbit polyclonal	WB, ELISA, IHC	affinity purified
FABP3 Human	Sheep polyclonal	ELISA	affinity purified
Leptin Human	Rabbit polyclonal	WB, ELISA, IP	affinity purified
NT-proANP Human	Sheep polyclonal	ELISA	affinity purified
Procalcitonin Human	Rabbit polyclonal	ELISA	affinity purified
	Sheep polyclonal	ELISA	affinity purified
Uromodulin Human	Rabbit polyclonal	WB, ELISA, IHC, IP	affinity purified
	Sheep polyclonal	WB, ELISA, IHC, IP	affinity purified
CTGF Human	Rabbit polyclonal	ELISA	affinity purified
	Sheep polyclonal	ELISA	affinity purified
PCSK9 Human	Rabbit polyclonal	ELISA	affinity purified
	Sheep polyclonal	ELISA	affinity purified
OSF-2 Human	Rabbit polyclonal	WB, IHC, ICC	affinity purified
HBx Virus Human	Rabbit polyclonal	WB, IHC, ICC, IP	affinity purified
Cathepsin Human	Mouse monoclonal, Clone: 3F9	WB, FC, IHC, IP	from serum free supernatant

* tested for named applications only

Date of issue October 2018

» Contact Information

BioVendor – Laboratorní medicína a.s.

Karasek 1767/1, 621 00 Brno, Czech Republic

+420 549 124 185 · +420 549 211 460 · E-mail: info@biovendor.com

www.biovendor.com

