

IMMUNO-TROL

A comprehensive range menu of internal and independent auto-immune Quality Control sera for FEIA-based systems



MULTI parameters QC for the FEIA-based systems



For assay validation

Universal multi-parametric controls

A comprehensive test menu consolidation, up to 6 parameters per vial. Reduces Quality Control storage, is cost effective, saves time and improves efficiency.

Internal and Independent quality controls

An unbiased performance assesment with FEIA-based system, Third Party QC complying with the requirements of ISO 15189.

Commutability

Identical use and reactivity as patient sample fulfilling the requirements of ISO 15189 and reducing the cost in Quality Control when reagent batch is changed.

Standardised controls

Our manufacturing processes guarantee the lot won't change frequently, the parameters levels will be mostly preserved lot-to-lot.

Shelf life


Shelf life up to two years. The same QC batch for a long time with limited need of QC target adjustments.

Liquid control

Ready-to-use. More accuracy and precision.

Service & Support

Do you have any question? Contact our Customer Service.



A large Quality Control range in Auto-immunity

IMMUNO-TROL  28 parameters

	dsDNA	U1-RNP	SmDp	Ro52	Ro60 (SS-A)	La (SS-B)	CENP	ScF-70S	Jo-1	Rib-P	Cardiolipin IgG	β2 GPI IgG	Cardiolipin IgM	β2 GPI IgM	CCP	RF IgM	TPO	TG	TSH-R*	MPO	PR3	GBM	Celikey IgG	Gliadin DP IgG	Celikey IgA	Gliadin DP IgA	ASCA IgA	ASCA IgG	
CONNECTIVE 1 FEIA	■	■	■	■	■	■																							
CONNECTIVE 2 FEIA							■	■	■	■																			
APS IgG FEIA											■	■																	
APS IgM FEIA												■	■																
RHEUMA FEIA														■	■														
THYRO FEIA																■	■	■											
VASCULITIS FEIA																			■	■	■								
CELIAC IgG FEIA																						■	■						
CELIAC IgA FEIA																								■	■				
ASCA FEIA																										■	■		

*Qualitative result, target to be established by each laboratory.

Product Description	Reference	Packaging
Connective 1 FEIA	HM 260-5X1	5X1 mL
Connective 2 FEIA	HM 261-5X1	5X1 mL
APS IgG FEIA	HM 262-5X1	5X1 mL
APS IgM FEIA	HM 263-5X1	5X1 mL
Rheuma FEIA	HM 264-5X1	5X1 mL
Thyro FEIA	HM-265-5X1	5X1 mL
Vasculitis FEIA	HM 267-5X1	5X1 mL
Celiac IgG FEIA	HM 269-5X1	5X1 mL
Celiac IgA FEIA	HM 270-5X1	5X1 mL
ASCA FEIA	HM 271-5X1	5X1 mL

FEATURES

- Human based matrix
- Liquid ready-to-use
- 28 parameters for 10 references
- Multi-parametric controls

SPECIFICATIONS

- 60 days open vial stability at +2°C to +8°C
- 12 months closed vial stability at +2°C to +8°C for IgM and 24 months for IgG and IgA

INSTRUMENTS

- IMMUNO-TROL to be used on FEIA-based systems

