

RENAL DISEASE

IMMUNOASSAYS

Cat. No.	Description	Assay Format	Sample Type	Assay Range
RD195023100	Adiponectin Human ELISA (Competitive)	Competitive ELISA, Immobilized antigen	Serum, Plasma	0.1 - 10 µg/ml
RD191023100	Adiponectin Human ELISA, High Sensitivity (Sandwich)	Sandwich ELISA, HRP-labelled antibody	Serum, Urine, Cerebrospinal fluid, Milk, Plasma	
RD293023100R	Adiponectin Mouse ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum	0.25 - 8 ng/ml
RSHAKMAN011R	Adiponectin HMW Mouse/Rat ELISA	Sandwich ELISA, HRP-labelled antibody	Serum	3.13 - 200 ng/ml
RSHAKRBS018R	Albumin Bovine ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma (EDTA, heparin, citrate), Cell culture supernatant	0.78 - 50 ng/ml
RSHAKRAL121R	Albumin Mouse ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Cell culture medium, Urine, Plasma	50 - 1000 ng/ml
RSHAKRAL120R	Albumin Rat ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma (EDTA, heparin, citrate), Cell culture supernatant, Urine	50 - 1000 ng/ml
RA19001R	Albumin Rat ELISA High Sensitivity	Competitive ELISA	Urine	8 - 1000 ng/ml
RSHAKRAL020SR	Albumin Urinary (S-type) Rat Assay	Immunturbidimetry/nephelometry	Serum, Urine, Plasma	6.17 - 500 µg/ml
RSHAKRAL021SR	Albumin Urinary Mouse (S-type) Assay	Immunturbidimetry/nephelometry	Serum, Urine	6.17 - 500 µg/ml
RA05880R	Angiotensin II Human/Rat ELISA	Sandwich ELISA, AChE-labelled antibody	Plasma (EDTA)	0.98 - 125 pg/ml
RD191367100R	Angiotensinogen Human ELISA	Sandwich ELISA	Serum, Urine, Cerebrospinal fluid, Bronchoalveolar lavage, Amniotic fluid, Plasma	0.31 - 20 ng/ml

RENAL DISEASE

Cat. No.	Description	Assay Format	Sample Type	Assay Range
RD191371200R	Annexin A1 Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Urine, Saliva, Plasma	0.1 – 3.2 ng/ml
RCD003R	Beta2-Microglobulin Human ELISA	Competitive ELISA, Immobilized antigen	Serum	0.2 – 10 mg/l
RSCPRH111R	Beta2-Microglobulin Rat ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma (EDTA, heparin, citrate), Urine	0.31 – 10 ng/ml
RD491034200R	Clusterin Canine ELISA (ApoJ)	Sandwich ELISA, Biotin-labelled antibody	Serum, Urine	4 – 128 ng/ml
RD194034200R	Clusterin Human ELISA (ApoJ)	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma, Urine, CSF, Cell lysate	5 – 160 ng/ml
RD191035200R	Connective Tissue Growth Factor Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (EDTA, citrate, heparin), Urine	0.63 – 20 ng/ml
RD491009100R	Cystatin C Canine ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Urine	0.31 – 10 ng/ml
RD191009100	Cystatin C Human ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma, Urine, CSF	200 – 10000 ng/ml
RD291009200R	Cystatin C Mouse ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum	0.156 – 10 ng/ml
RD391009200R	Cystatin C Rat ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum	0.078 – 25 ng/ml
RAF013R	Erythropoietin Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (EDTA, heparin, citrate), Cell culture supernatant	1.6 – 100 mIU/ml
RD191037100	Fetuin-A Human ELISA (AHSN)	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma	2 – 100 ng/ml
RD191407100R	Haptoglobin Human ELISA	Sandwich ELISA, HRP-labelled antibody	Cerebrospinal fluid, Plasma (citrate, EDTA, heparin), Serum, Urine	3.13 – 200 ng/ml
RD192247200R	Heart FABP (FABP3) Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Urine, Cerebrospinal fluid, Bronchoalveolar lavage, Amniotic fluid, Plasma	20 – 1280 pg/ml
RD191422200R	IGFBP-7 Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Tissue homogenates, Serum, Urine, Cerebrospinal fluid, Amniotic fluid, Plasma	0.63 – 40 ng/ml
RD191102200R	Lipocalin 2/ NGAL Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (EDTA), Urine, Stool	0.3 – 10 ng/ml
RD191486200R	NT-proBNP Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Urine, Cerebrospinal fluid, Plasma	5 – 230 pg/ml
RD191446200R	Osteopontin Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Cerebrospinal fluid, Milk, Serum, Plasma, Urine	0.25 – 16 ng/ml
RD191390200R	Placental Growth Factor Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Urine, Cerebrospinal fluid, Amniotic fluid, Plasma	2 – 128 pg/ml
RMI001R	Prorenin Human ELISA	Direct ELISA, HRP-labelled antibody	Plasma	0 – 10 ng/ml
RD191167200R	Total Prorenin/Renin Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Amniotic fluid, Plasma (EDTA, citrate, heparin) Serum, Urine	25 – 1600 pg/ml
RD191113100R	Prostaglandin D Synthase (Lipocalin-type) Human ELISA	Sandwich ELISA, HRP-labelled antibody	Serum, Plasma, Urine, CSF	1 – 40 ng/ml
RD191069200R	Prouroguanylin Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma (EDTA, heparin, citrate)	0.8 – 40 ng/ml
RD194169200R	Active Renin Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Amniotic fluid, Plasma (EDTA, citrate, heparin), Serum, Urine	4.5 – 144 μ U/ml

RENAL DISEASE

Cat. No.	Description	Assay Format	Sample Type	Assay Range
RDS004R	Renin Human ELISA	Sandwich ELISA, HRP-labelled antibody	Plasma (EDTA)	0 - 270 pg/ml
RD191116200R	sRAGE Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma	50 - 3200 pg/ml
RD191408200R	suPAR Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma, Urine	7.8 - 500 pg/ml
RD491163200R	Uromodulin Canine ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Urine	0.16 - 10 ng/ml
RD191163200R	Uromodulin Human ELISA	Sandwich ELISA, Biotin-labelled antibody	Serum, Plasma, Urine	0.5 - 32 ng/ml

RENAL DISEASE

PROTEINS

Cat. No.	Description	Species	Type	Source
RD162023050	Adiponectin	Human	Native	Human pooled serum
RD172023025-F	Adiponectin	Human	Recombinant	HEK293
RD172023100	Adiponectin	Human	Recombinant	HEK293
RD172029100	Adiponectin	Human	Recombinant	<i>E. coli</i>
RD272023100	Adiponectin	Mouse	Recombinant	HEK293
RD572023100	Adiponectin	Porcine	Recombinant	HEK293
RD372023100	Adiponectin	Rat	Recombinant	HEK293
RD172112100	Adiponectin Globular	Human	Recombinant	<i>E. coli</i>
RD172023100-C	Adiponectin HMW-rich	Human	Recombinant	HEK293
RD172023100-F	Adiponectin Solution	Human	Recombinant	HEK293
RD272587025	Adiponectin, Tag-free	Mouse	Recombinant	HEK293
RD272587100	Adiponectin, Tag-free	Mouse	Recombinant	HEK293
RD172091025	Adiponectin, Trimeric form (HEK)	Human	Recombinant	HEK293
RD172091100	Adiponectin, Trimeric form (HEK)	Human	Recombinant	HEK293
RD272091100	Adiponectin, Trimeric form (HEK)	Mouse	Recombinant	HEK293
RD172367100	Angiotensinogen	Human	Recombinant	HEK293
RD172371100	Annexin A1	Human	Recombinant	<i>E. coli</i>
RD162034025	Clusterin	Human	Native	Serum
RD472034100	Clusterin (Apo J)	Canine	Recombinant	<i>E. coli</i>
RD472034100-HEK	Clusterin (Apo J)	Canine	Recombinant	HEK293
RD172034100	Clusterin (Apo J)	Human	Recombinant	HEK293
RD172034050-S	Clusterin (Apo J)	Human	Native	Plasma
RD372034100	Clusterin (Apo J)	Rat	Recombinant	<i>E. coli</i>
RD176034100	Clusterin, Tag free	Human	Recombinant	HEK293
RD172035100	Connective Tissue Growth Factor	Human	Recombinant	<i>E. coli</i>
RD172035100-HEK	Connective Tissue Growth Factor	Human	Recombinant	HEK293
RD272589100	Connective Tissue Growth Factor	Mouse	Recombinant	HEK293
RD176035025	Connective Tissue Growth Factor, Tag Free	Human	Recombinant	HEK293

RENAL DISEASE

Cat. No.	Description	Species	Type	Source
RD172591100	CTGF C-terminus	Human	Recombinant	<i>E. coli</i>
RD172593100	CTGF, TSP Type-1 Domain	Human	Recombinant	<i>E. coli</i>
RD472009100	Cystatin C	Canine	Recombinant	<i>E. coli</i>
RD172009100	Cystatin C	Human	Native	Urine
RD172009100-H	Cystatin C	Human	Recombinant	<i>E. coli</i>
RD272009100	Cystatin C	Mouse	Recombinant	<i>E. coli</i>
RD372009100	Cystatin C	Rat	Recombinant	<i>E. coli</i>
RD172037100	Fetuin-A (AHSG)	Human	Native	Plasma
RD172241100	Fibroblast Growth Factor 23 C-terminal peptide	Human	Recombinant	<i>E. coli</i>
RD172150100	KIM-1	Human	Recombinant	HEK293
RD172099100	Klotho	Human	Recombinant	CHO
RD172370100	Leukocyte Cell-Derived Chemotaxin-2	Human	Recombinant	<i>E. coli</i>
RD172102050	Lipocalin-2/NGAL	Human	Recombinant	<i>E. coli</i>
RD172240100	Liver FABP	Human	Recombinant	<i>E. coli</i>
RD172240101	Liver FABP, Tagless	Human	Recombinant	<i>E. coli</i>
RD172170100	Osteoactivin/GPNMB	Human	Recombinant	HEK293
RD172446100	Osteopontin	Human	Recombinant	HEK293
RD172330100	Podocin	Human	Recombinant	<i>E. coli</i>
RD172046100	Proguanylin	Human	Recombinant	<i>E. coli</i>
RD172113100	Prostaglandin-D Synthase (Lipocalin-type)	Human	Recombinant	<i>E. coli</i>
RD172069100	Prouroguanylin	Human	Recombinant	<i>E. coli</i>
RD172116100	sRAGE	Human	Recombinant	<i>E. coli</i>
RD172116100-HEK	sRAGE	Human	Recombinant	HEK293
RD172167100-HEK	Renin	Human	Recombinant	HEK293
RD172264100	Renin Receptor	Human	Recombinant	<i>E. coli</i>
RD172116025-HEK	Soluble Advanced Glycation End Product-Specific Receptor	Human	Recombinant	HEK293
RD172160100	Trefoil Factor 3	Human	Recombinant	<i>E. coli</i>
RD372160100	Trefoil Factor 3	Rat	Recombinant	<i>E. coli</i>
RD463163100	Uromodulin	Canine	Native	Canine urine
RD663163050	Uromodulin	Feline	Native	Feline urine

RENAL DISEASE

Cat. No.	Description	Species	Type	Source
RD172163100	Uromodulin	Human	Native	Human urine
RD563163100	Uromodulin	Porcine	Native	Porcine urine

ANTIBODIES

Cat. No.	Description	Spec. React.	Type	Antigen
RD1820231005E4	Adiponectin	Human	Mouse MAb, Clone: 5E4	HEK293
RD1820231007C11	Adiponectin	Human	Mouse MAb, Clone: 7C11	HEK293
RD181023220	Adiponectin	Human	Rabbit PAb	<i>E. coli</i>
RD184023100	Adiponectin (HEK)	Human	Goat PAb	HEK293
RD181023100	Adiponectin (HEK)	Human	Rabbit PAb	HEK293
RD184023100-01	Adiponectin (HEK)	Human	Sheep PAb	HEK293
RD581023100	Adiponectin (HEK)	Porcine	Rabbit PAb	HEK293
RD584023100	Adiponectin (HEK)	Porcine	Sheep PAb	HEK293
RD181371100	Annexin A1	Human	Rabbit PAb	<i>E. coli</i>
RD184371100	Annexin A1	Human	Sheep PAb	<i>E. coli</i>
RD184328100	Apoptosis Inhibitor of Macrophage	Human	Sheep PAb	<i>E. coli</i>
RD481034100	Clusterin (Apo J)	Canine	Rabbit PAb	<i>E. coli</i>
RD484034100	Clusterin (Apo J)	Canine	Sheep PAb	<i>E. coli</i>
RD182034110-H3	Clusterin (Apo J)	Human	Mouse MAb, Clone Hs-3	Human sperm extract
RD1820341008H4-01	Clusterin	Human	Mouse MAb, Clone: 8H4	Human plasma
RD381034100	Clusterin (Apo J)	Rat	Rabbit PAb	<i>E. coli</i>
RD181035100	Connective Tissue Growth Factor	Human	Rabbit PAb	HEK293
RD184035100	Connective Tissue Growth Factor	Human	Sheep PAb	HEK293
RD481009100	Cystatin C	Canine	Rabbit PAb	<i>E. coli</i>
RD484009100	Cystatin C	Canine	Sheep PAb	<i>E. coli</i>
RD181009220-01	Cystatin C	Human	Rabbit PAb	Human urine
RD181009100	Cystatin C	Human	Rabbit PAb	<i>E. coli</i>
RD184009100	Cystatin C	Human	Sheep PAb	<i>E. coli</i>

RENAL DISEASE

Cat. No.	Description	Spec. React.	Type	Antigen
RD281009100	Cystatin C	Mouse	Rabbit PAb	<i>E. coli</i>
RD284009100	Cystatin C	Mouse	Sheep PAb	<i>E. coli</i>
RD384009100	Cystatin C	Rat	Sheep PAb	<i>E. coli</i>
RD184036100	Fetuin-A (AHSG)	Human	Goat PAb	Human serum
RD181036100	Fetuin-A (AHSG)	Human	Rabbit PAb	Human serum
RD181109100	Fibroblast Growth Factor 23	Human	Rabbit PAb	<i>E. coli</i>
RD184240100	Liver Fatty Acid Binding Protein	Human	Sheep PAb	<i>E. coli</i>
RD181170100	Osteoactivin/GPNMB	Human	Rabbit PAb	HEK293
RD184170100	Osteoactivin/GPNMB	Human	Sheep PAb	HEK293
RD181113100	Prostaglandin-D Synthase (Lipocalin-type)	Human	Rabbit PAb	<i>E. coli</i>
RD184113100	Prostaglandin-D Synthase (Lipocalin-type)	Human	Sheep PAb	<i>E. coli</i>
RD184116100	sRAGE	Human	Sheep PAb	<i>E. coli</i>
RD181221100	S100A12	Human	Rabbit PAb	<i>E. coli</i>
RD184221100	S100A12	Human	Sheep PAb	<i>E. coli</i>
RD181160100	Trefoil Factor 3	Human	Rabbit PAb	<i>E. coli</i>
RD184160100	Trefoil Factor 3	Human	Sheep PAb	<i>E. coli</i>
RD181343100	Urinary Trypsin Inhibitor	Human	Rabbit PAb	Human Urine
RD181163100	Uromodulin	Human	Rabbit PAb	<i>E. coli</i>
RD184163100	Uromodulin	Human	Sheep PAb	<i>E. coli</i>

› Contact Information



BioVendor – Laboratorni medicina a.s.
Karasek 1767/1, 621 00 Brno, Czech Republic
Phone: +420 549 124 185
Fax: +420 549 211 460
E-mail: info@biovendor.com

BioVendor GmbH
Otto-Hahn-Straße 16, 34123 Kassel, Germany
Phone: +49 (0)561 861 509 33
Fax: +49 (0)561 861 509 34
E-mail: infoEU@biovendor.com

› www.biovendor.com

BioVendor GesmbH
Gaudenzdorfer Gürtel 43-45, 1120 Vienna, Austria
Phone: +43 1 890 9025
Fax: +43 1 890 5163
E-mail: infoAustria@biovendor.com

BioVendor, LLC
128 Bingham Rd., Suite 1300, Asheville, NC 28806, USA
Phone: +1-800-404-7807, Phone: +1-828-575-9250
Fax: +1-828-575-9251
E-mail: infoUSA@biovendor.com



Oxford Biosystems Ltd
184B Park House, Park Drive
Milton Park, Abingdon, Oxfordshire
OX14 4SR, United Kingdom
Phone: +44 (0)1235 431390
E-mail: sales@oxfordbiosystems.com

› www.oxfordbiosystems.com